



C1

reusable fractional electrodes catalog




reusable fractional electrodes



64 Pins	Wavetronic 6000 Touch	Wavetronic 5000 Digital
Reference	ACELO907	ACELO255

Square fractional tip, 64 pins, 8 x 8
Ø 0,2 mm



36 Pins	Wavetronic 6000 Touch	Wavetronic 5000 Digital
Reference	ACELO911	ACELO257

Square fractional tip, 36 pins, 6 x 6
Ø 0,2 mm




16 Pins	Wavetronic 6000 Touch	Wavetronic 5000 Digital
Reference	ACELO915	ACELO259

Square fractional tip, 16 pins, 4 x 4
Ø 0,2 mm

Description

Electrode with fractional tip for megapulse HF Fraxx, autoclave compatible, cable and connector included

The fractional electrode when contact the skin produce micro lesions on dermis without cause damage on, because its not invasive, the fractional energy will be distributed, keeping the adjacent tissues untouched. This process will stimulate the production of collagen and the regeneration of dermis on the interface level between papillary and reticular dermis through fibroblast stimulation.




24 Pins	Wavetronic 6000 Touch	Wavetronic 5000 Digital
Reference	ACELO925	ACELO258

Rectangular fractional tip, 24 pins, 8 x 3
Ø 0,2 mm

Treatments

- Flaccid Skin
- Wrinkles
- Dark Circles
- Skin Spots
- Acne Scar
- Hypertrophy Scar
- Facial Rejuvenation
- Body Rejuvenation
- Stretch mark
- And more



24 Pins	Wavetronic 6000 Touch	Wavetronic 5000 Digital
Reference	ACELO931	ACELO263

Rectangular fractional tip, 24 pins, 8 x 3
Ø 0,2 mm